



MEDICA DEPOT

**JUVÉDERM® ULTRA
PLUS XC**
Product Specifications

JUVÉDERM® ULTRA PLUS XC

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed physician or properly licensed practitioner.

BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY.

1. DEVICE DESCRIPTION

JUVÉDERM® Ultra Plus XC is a sterile, biodegradable, non-pyrogenic, viscoelastic, clear, colorless, homogeneous gel implant. It consists of cross-linked hyaluronic acid (HA) produced by *Streptococcus* species of bacteria, formulated to a concentration of 24 mg/mL and 0.3% w/w lidocaine in a physiologic buffer.

2. INTENDED USE/INDICATIONS

JUVÉDERM® Ultra Plus XC injectable gel is indicated for injection into the mid to deep dermis for correction of moderate to severe facial wrinkles and folds (such as nasolabial folds).

3. CONTRAINDICATIONS

- JUVÉDERM® Ultra Plus XC is contraindicated for patients with severe allergies manifested by a history of anaphylaxis or history or presence of multiple severe allergies.
- JUVÉDERM® Ultra Plus XC contains trace amounts of Gram-positive bacterial proteins and is contraindicated for patients with a history of allergies to such material.
- JUVÉDERM® Ultra Plus XC contains lidocaine and is contraindicated for patients with a history of allergies to such material.

4. WARNINGS

- The product must not be injected into blood vessels. Introduction of JUVÉDERM® Ultra Plus XC into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Take extra care when injecting soft tissue fillers; for example, after insertion of the needle, and just before injection, the plunger rod can be withdrawn slightly to aspirate and verify the needle is not intravascular, inject the product slowly and apply the least amount of pressure necessary. Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage leading to stroke, skin necrosis, and damage to underlying facial structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur (see Health Care Professional Instructions #11).
- Product use at specific sites in which an active inflammatory process (skin eruptions such as cysts, pimples, rashes, or hives) or infection is present should be deferred until the underlying process has been controlled.
- Injection procedure reactions consist mainly of short-term inflammatory symptoms starting early after treatment and lasting ≤ 7 days' duration. Refer to the ADVERSE EVENTS section for details.

5. PRECAUTIONS

- JUVÉDERM® Ultra Plus XC is packaged for single-patient use. Do not resterilize. Do not use if package is opened or damaged.
- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.

- Health care professionals are encouraged to discuss all potential risks of soft-tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Based on preclinical studies, patients should be limited to 20 mL of any JUVÉDERM® injectable gel per 60 kg (132 lbs) body mass per year. The safety of injecting greater amounts has not been established.
- The safety and effectiveness for the treatment of anatomic regions other than facial wrinkles and folds (eg, lips) have not been established in controlled clinical studies.
- As with all transcutaneous procedures, dermal filler implantation carries a risk of infection. Standard precautions associated with injectable materials should be followed.
- JUVÉDERM® Ultra Plus XC is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity and performance of the product and it can therefore no longer be assured.
- The safety for use during pregnancy, in breastfeeding females, or in patients under 18 years has not been established.
- The safety in patients with known susceptibility to keloid formation, hypertrophic scarring and pigmentation disorders has not been studied.
- JUVÉDERM® Ultra Plus XC should be used with caution in patients on immunosuppressive therapy.
- Patients who are using substances that can prolong bleeding (such as aspirin, nonsteroidal anti-inflammatory drugs, and warfarin) may, as with any injection, experience increased bruising or bleeding at injection sites.
- After use, treatment syringes and needles may be potential biohazards. Handle and dispose of these items in accordance with accepted medical practice and applicable local, state, and federal requirements.
- JUVÉDERM® Ultra Plus XC injectable gel is a clear, colorless gel without visible particulates. In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe; notify Allergan® Product Support at 1-877-345-5372.
- If laser treatment, chemical peeling or any other procedure based on active dermal response is considered after treatment with JUVÉDERM® Ultra Plus XC, there is a possible risk of eliciting an inflammatory reaction at the implant site. An inflammatory reaction is also possible if the product is administered before the skin has healed completely after such a procedure.
- Failure to comply with the needle attachment instructions could result in needle disengagement and/or product leakage at the LUER-LOK® and needle hub connection.

6. ADVERSE EVENTS

A. Clinical Evaluation of JUVÉDERM® Ultra Plus XC

A 2-week, randomized, controlled US clinical study for JUVÉDERM® Ultra XC and Ultra Plus XC compared with JUVÉDERM® Ultra and Ultra Plus without lidocaine showed a similar safety profile in all subjects (N = 72), with the exception of fewer reports of pain/tenderness with the product containing lidocaine. Common treatment site responses (CTR) by severity and duration, are presented in Tables 1 and 2. Aside from injection site responses, there were no adverse events related to the device, procedure, or anesthesia.

- The most common injection site responses for JUVÉDERM® Ultra Plus XC were redness, swelling, tenderness, firmness, lumps/bumps, discoloration, and bruising.

Table 1. Injection Site Responses by Maximum Severity (Number/% of Subject Nasolabial Folds [NLFs])

Injection Site Responses	JUVEDERM® Ultra Plus XC (N = 38 NLFs)				JUVEDERM® Ultra Plus (N = 38 NLFs)			
	JUVEDERM® Ultra Plus XC n	JUVEDERM® Ultra Plus XC n	Mild n	Severe n	JUVEDERM® Ultra Plus n	JUVEDERM® Ultra Plus n	Mild n	Severe n
Redness	27	26	19	7	1	18	8	0
Pain	17	24	11	6	0	15	5	4
Tenderness	29	43	21	11	3	42	11	13
Firmness	31	34	21	7	3	21	10	3
Swelling	86%	94%	59%	19%	6%	58%	28%	6%
Lumps/Bumps	27	27	19	9	1	15	11	1
Itching	28	28	18	6	3	42%	31%	3%
Discoloration	28	28	18	6	1	11%	44%	31%
	28	28	20	5	3	31%	0%	5%
	28	28	20	5	3	31%	0%	5%
	28	28	20	5	3	31%	0%	5%

^a Number of subject NLFs treated with the respective device
^b Mild = Moderate
^c Number of NLFs with any occurrence of a particular CTR (or severity for the overall percentage)

Table 3. Injection Site Responses by Maximum Severity (Number/% of Subject NLFs)

Injection Site Responses	TOTALS				JUVEDERM® Ultra Plus (N = 144 NLFs)				ZYPLAST™ Ultra Plus (N = 144 NLFs)			
	JUVEDERM® Ultra Plus n	JUVEDERM® Ultra Plus n	Mild n	Severe n	JUVEDERM® Ultra Plus n	JUVEDERM® Ultra Plus n	Mild n	Severe n	ZYPLAST™ Ultra Plus n	ZYPLAST™ Ultra Plus n	Mild n	Severe n
Redness	129	128	61	61	7	42	15	10%	10%	42	15	10%
Pain	97	89	42%	46%	5	38	22	5	38	22	5	38
Tenderness	127	122	59	53	10	61	22	3	61	22	3	61
Firmness	88%	85%	41%	37%	10%	43%	35%	6%	43%	35%	6%	43%
Swelling	124	121	61	50	13	71	41	9	71	41	9	71
Lumps/Bumps	70	68	42%	38%	9%	48%	28%	8%	48%	28%	8%	48%
Itching	48	49	51	38	9	2	39	9	39	2	39	9
Discoloration	34%	33%	26%	20%	10%	3%	22%	6%	26%	10%	3%	22%

^a Number of subject NLFs treated with the respective device
^b Mild = Moderate
^c Number of subject NLFs with each specific injection site response

unsatisfactory result, allergic reaction, skin discoloration, vascular occlusion, device migration, infection, neurological symptoms such as increase or decrease in sensation, inflammatory nodule/granuloma, dermatitis, blister, dry skin, anxiety, overcorrection, necrosis, bleeding, abscess, herpes, flu-like symptoms, scarring, varied injection, angioedema, vision abnormalities, acne, headache, malaise, drainage, dyspnea, extrusion, cyst, dizziness, syncope, telangiectasia, anaphylactic reaction, calcification, depression, nausea, autoimmune disorder exacerbation, beading, deeper wrinkle, cardiac complications, vision loss, and traumatic injury.

In many cases, the symptoms resolved without any treatment. Reported treatments have included: antibiotics, steroids, steroid creams, hyaluronidase, anti-inflammatories, anti-histamines, need aspiration and drainage, hyperbaric oxygen, laser resurfacing, anti-viral, excision, eye drops, hyperbaric oxygen, laser resurfacing, tissue debridement, surgical scar revision, ice, massage, jelly, anxiolytics, antifungals, anticoagulants, and epinephrine.

Inflammatory reaction at the injection site, mostly a nonserious event, has been reported in association with edema, erythema, ecchymosis, pruritus, induration, pain, nodule, blister, abscess, and infection. Time to onset ranged from 1 day to 4 months post JUVEDERM® Ultra Plus injection, and outcomes ranged from the health care professionals included topical steroid cream, oral steroids, and antibiotics. Additional treatment noted was a needle aspiration for drainage of an abscess.

Vascular occlusion of vessels resulting in necrosis and vision abnormalities have been reported following injection of JUVEDERM® products, with and without lidocaine, with a time to onset ranging from immediate to within 1 week following injection. These reported events likely resulted from inadvertent arterial injection. In many of these cases, the product was injected into the highly vascularized areas of the glabella, nose, and periorbital area, which are outside the device indications for use (see WARNINGS section). Reported treatments include: anticonvulsants, epinephrine, aspirin, hyaluronidase, steroid treatment, eye drops, hyperbaric oxygen, surgery, vasodilators, and warm compress. Outcomes have ranged from completely resolved to ongoing at the time of last contact.

Serious adverse events have infrequently been reported for JUVEDERM® Ultra Plus (reported with a frequency of 5 or more). The most commonly reported serious adverse events were edema, erythema, ecchymosis, and pain.

- The onset of edema, erythema, and pain generally varied from immediate to 2 months post injection. The treatment prescribed included NSAIDs, anti-histamines, antibiotics, steroids, and hyaluronidase. In most cases, the reported events resolved within a few days to 5 weeks.
- The onset of ecchymosis generally varied from immediate to 1 day post injection. The treatment prescribed included NSAIDs, anti-histamines, antibiotics, steroids, and hyaluronidase. In most cases ecchymosis resolved within a few days to 6 weeks.

Additionally, there have been reports of nodules, infection, and inflammation.

- The onset of nodules generally varied from immediate to 2 months post injection. The treatment prescribed included NSAIDs, antibiotics, steroids, and hyaluronidase. In most cases nodules resolved within 1 month.
- The onset of infection generally varied from immediate to 1 month post injection. The treatment prescribed included antibiotics, pain killers, and antibacterial drugs.

The onset of inflammation generally varied from the day of treatment to 1 day post injection. The treatment prescribed included antibiotics, steroids, and needle aspiration. Resolution of symptoms has been reported within 4 days.

Delayed-onset inflammation near the site of dermal filler injections is one of the known adverse events associated with dermal fillers. Cases of delayed-onset inflammation have been reported to occur at the dermal filler treatment site following viral or bacterial illnesses or infections, vaccinations, or dental procedures. Typically, the reported inflammation was responsive to treatment or resolved on its own.

Adverse reactions should be reported to the Allergan® Product Surveillance Department at 1-877-345-5372.

7. CLINICAL STUDIES

A. Pivotal Study for JUVEDERM® Ultra Plus (Without Lidocaine)

A prospective, double-blind, randomized, within-subject, controlled, multicenter pivotal clinical study was conducted to evaluate the safety and effectiveness of JUVEDERM® Ultra Plus in the treatment of moderate to severe wrinkles. Subjects underwent treatment with JUVEDERM® Ultra Plus in one NLF and the control implant (ZYPLAST™ bovine collagen) in the opposite NLF.

Up to 3 bilateral treatments (initial treatment and up to 2 touch-up treatments), approximately 2 weeks apart, were allowed. At 2 and 4 weeks after each treatment, the independent Expert Reviewer (IER) assessed the level of correction achieved. If correction was less than optimal after the first or second treatment, the investigator re-treated the under-corrected NLFs using the same respective treatment materials as in the initial treatment. The IER and the subject remained masked to the randomized treatment assignment.

Routine follow-up visits for safety and effectiveness occurred at days 3 and 7 and week 2 after each treatment, and at 4, 8, 12, 16, 20, and 24 weeks after the last treatment. Standardized facial photography was performed for documentation purposes. The investigator and the IER independently evaluated the severity of the subject's NLFs using a validated 5-point (range 0 to 4) photographic NLF severity scale. The subject made independent self-assessments of NLF severity using a non-photographic 5-point grading scale.

Study Endpoints

The primary effectiveness endpoint for the study was the IER's NLF severity score over the post-treatment follow-up period. Effectiveness of device treatment was demonstrated by a lowering of the NLF severity score. Additional analyses included the subject's and the investigator's live NLF severity assessments.

Subject Demographics

A total of 146 subjects (26 to 74 years of age) were randomized and treated, and 140 (96%) completed the 6-month follow-up period. Prior to enrollment, 80 (55%) had previous experience with other facial dermal treatments (eg, alpha-hydroxy agents, neurotoxins, microdermabrasion or retinoic acid).

Subject demographics and pretreatment characteristics of the JUVEDERM® Ultra Plus effectiveness population are presented in Table 5.

Table 4. Duration of Injection Site Responses Occurring in > 5% of Subject NLFs

Injection Site Responses	TOTALS				JUVEDERM® Ultra Plus (N = 144 NLFs)				ZYPLAST™ Ultra Plus (N = 144 NLFs)			
	JUVEDERM® Ultra Plus n	JUVEDERM® Ultra Plus n	Mild n	Severe n	JUVEDERM® Ultra Plus n	JUVEDERM® Ultra Plus n	Mild n	Severe n	ZYPLAST™ Ultra Plus n	ZYPLAST™ Ultra Plus n	Mild n	Severe n
Duration ^a	≤ 3 Days	4-7 Days	8-14 Days	≥ 3 Days	4-7 Days	8-14 Days	≥ 3 Days	4-7 Days	8-14 Days	≥ 3 Days	4-7 Days	≥ 14 Days
Redness	10	30%	7%	14%	37%	26%	9%	17%	10	30%	7%	14%
Pain	59	37	25	8	55	34	17	7	59	37	25	8
Tenderness	41%	26%	17%	6%	39%	31%	12%	5%	41%	26%	17%	6%
Firmness	24	28	18	56	28	26	16	52	24	28	18	56
Swelling	31	40	21	52	32	52	12	33	31	40	21	52
Lumps/Bumps	32	24	19	45	15	37%	33%	9%	6%	32	24	19
Itching	25	17%	13%	31%	10%	18%	10%	40%	25	17%	13%	31%
Discoloration	22	11	4	12	27	13%	6%	2%	22	11	4	12

^a Number of subject NLFs treated with the respective device
^b Number of subject NLFs with each specific injection site response by maximum duration of inflammation
^c Duration refers to number of days from symptom onset until resolution, irrespective of date of implantation

Local injection site responses were recorded in subjects' diaries one or more times for 98% of JUVEDERM® Ultra Plus treated subjects and 97% of ZYPLAST™ treated NLFs. Subjects' scores for both products were predominantly Mild or Moderate in intensity, and their injection site responses reported by greater than 1% of subjects and not noted in the above tables were skin tingling and peeling. No clinically meaningful differences in the safety profiles of JUVEDERM® Ultra Plus and ZYPLAST™ were found during the study.

C. Other Safety Data

In 2 additional randomized U.S. clinical studies of other JUVEDERM® formulations (without lidocaine) in a total of 295 subjects, the safety profile was similar to that described above for JUVEDERM® Ultra Plus.

Postmarket Surveillance

The following adverse events were received from postmarket surveillance for JUVEDERM Ultra and Ultra Plus, with and without lidocaine, which were not observed in the clinical trials; this includes reports received globally from all sources including scientific journals and voluntary reports. All adverse events obtained through postmarket surveillance with a frequency of 5 or more events are listed in order of prevalence: edema, inflammatory reaction, non-inflammatory nodule, lack of loss of correction, pain, hematoma,

Table 2. Duration of Injection Site Responses Occurring in > 5% of Subject NLFs

Injection Site Responses	TOTALS				JUVEDERM® Ultra Plus XC (N = 38 NLFs)				JUVEDERM® Ultra Plus (N = 38 NLFs)			
	JUVEDERM® Ultra Plus XC n	JUVEDERM® Ultra Plus XC n	Mild n	Severe n	JUVEDERM® Ultra Plus n	JUVEDERM® Ultra Plus n	Mild n	Severe n				
Duration ^a	1-3 Days	4-7 Days	≥ 14 Days	≥ 14 Days	1-3 Days	4-7 Days	≥ 14 Days	≥ 14 Days				
Redness	17	5	3	2	17	5	2	2				
Pain	16	14	0	0	24	0	0	0				
Tenderness	4%	3%	0%	67%	0%	0%	0%	0%				
Firmness	50%	22%	8%	53%	22%	8%	3%	3%				
Swelling	9	12	6	4	15	11	3	5				
Lumps/Bumps	25%	33%	17%	11%	42%	31%	8%	14%				
Itching	17	7	4	1	21	7	1	2				
Discoloration	4%	18%	1%	3%	56%	16%	3%	6%				
	28%	22%	6%	19%	33%	17%	3%	22%				
	9	7	9	3	12	11	3	2				
	25	19%	25%	8%	32%	31%	8%	6%				
	23	36	0	11	0	0	0	0				
	16	5	2	5	12	0%	0%	0%				
	4%	14%	6%	14%	33%	14%	6%	11%				

^a Number of subject NLFs treated with the respective device
^b Mild = Moderate
^c Duration refers to number of days from symptom onset until resolution, irrespective of date of implantation

B. Clinical Evaluation of JUVEDERM® Ultra Plus (Without Lidocaine)

In the initial randomized, controlled clinical trial to evaluate safety and effectiveness, 44 subjects were injected with JUVEDERM® Ultra Plus in one NLF and ZYPLAST™ dermal filler in the contralateral NLF. Preprinted diary forms were used by subjects to record specific signs and symptoms experienced during each of the first 14 days (day 0 through day 13) after initial and touch-up treatments. Subjects were instructed to rate each common treatment response listed on the diary as "Mild," "Moderate," "Severe," or "None." Injection site responses reported by > 5% of subjects in either treatment group are summarized in Tables 3 and 4.

Table 5. Demographics and Pretreatment Characteristics of the Effectiveness Population (Number/% of Subjects)
N = 146

Gender (Number/%)	Female	90%
Male	14	10%
Ethnicity (Number/%)	107	73%
Caucasian	17	12%
African American	20	14%
Hispanic	0	0%
Asian	0	0%
Other	2	1%
Fitzpatrick Skin Phototype (Number/%)	8	5%
I	34	23%
II	51	35%
III	31	21%
IV	18	12%
V	18	12%
VI	4	3%
Mean Baseline NLF Severity Score ^a	2.6	
JUVEDERM® Ultra Plus NLF	2.6	
ZYPLAST® NLF		

^a NLF Severity was ranked on a 5-point scale from None (0) to Extreme (4)

The primary effectiveness results for JUVEDERM® Ultra Plus based on the IER's assessment of NLF severity are presented in Table 6.

Table 6. Effectiveness Summary Independent Expert Reviewer's NLF Severity Scores

	JUVEDERM® Ultra Plus (N = 146 NLFs)	Control (N = 146 NLFs)
n ^a	146	146
NLF Severity ^b	2.6	2.6
Improvement Since Baseline ^c	-	-
Baseline	1.46	1.46
Week 2	0.9	0.9
Week 12	0.9	0.9
Week 24	1.2	1.4
Improvement Since Baseline ^d	2.2	0.4

^a Number of subject NLFs treated with the respective device

^b A commercially available legible Levine collagen implant

^c Number of subject NLFs with data at baseline and the specified time point

^d Mean score

Throughout the 24-week study period, JUVEDERM® Ultra Plus provided a clinically and statistically significant improvement in NLF severity. Clinical superiority was achieved at week 24 for JUVEDERM® Ultra Plus over ZYPLAST® with mean NLF severity of 1.2 and 2.2, respectively (P < 0.0001). Additionally, subject assessments for product preference overwhelmingly favored JUVEDERM® Ultra Plus; 84% preferred the JUVEDERM® Ultra Plus treated NLF over the ZYPLAST® treated NLF.

B. Extended Follow-up Clinical Study

Of the 146 randomized and treated subjects, more than three-quarters (76%, 111/146) returned after completion of their 24-week follow-up in the pivotal study for complimentary repeat treatment. Demographics for the subjects receiving repeat treatment were similar to those in the overall study. The majority of subjects were Caucasian and female, with a median age of 47 years. More than one-third of subjects were of Fitzpatrick Skin Phototypes IV, V or VI.

After completing the 24-week study, subjects returned for repeat treatment at their convenience or their investigator's convenience. The average time elapsed between last initial treatment and repeat treatment was approximately 9 months. A statistical analysis demonstrated that those subjects who returned for repeat treatment at a later time point were representative of the pivotal study subjects overall. There were no significant differences between these stratified groups in terms of NLF severity at baseline or at the 24-week follow-up visit or in overall initial volume injected. Before repeat treatment, live assessments of wrinkles severity were made by the investigator and the subject. The extended follow-up effectiveness results for JUVEDERM® Ultra Plus based on the investigator's assessment of NLF severity are presented in Table 7.

Table 7. Extended Follow-up Prior to Repeat Treatment Effectiveness Summary

	JUVEDERM® Ultra Plus (N = 111 NLFs)	Improvement Since Baseline ^a	P-value
Baseline ^b	111	2.6	N/A
Follow-up Week 24 ^c (Month 9)	110	1.1	< .0001
Follow-up Weeks 25-36 (Months 6-9)	64	1.1	< .0001
Follow-up Weeks 37-48 (Months 9-12)	24	1.4	< .0001
Follow-up Weeks >48 (> 1 year)	23	1.6	< .0001

^a Mean score

^b NLFs included during pivotal study

^c Mean score

All subjects returning for repeat treatment were stratified into 3 groups based on the time elapsed between last initial treatment and repeat treatment: 25 to 36 weeks, 37 to 48 weeks, or > 48 weeks. Mean improvement since baseline was clinically significant (> 1 point) for all 3 groups, with a large majority of subjects treated with JUVEDERM® Ultra Plus demonstrating improvement:

- 92% (59/64) at 25 to 36 weeks (6-9 months)
- 83% (20/24) at 37 to 48 weeks (9-12 months)
- 78% (18/23) beyond 48 weeks (beyond 1 year)

Follow-up After Repeat Treatment

A subset of subjects enrolled in a prospective, multicenter study for follow-up after repeat treatment. Subjects were eligible for the follow-up study if they completed the pivotal study, indicated that they preferred JUVEDERM® Ultra Plus over the control device, and received repeat treatment between 24 and 36 weeks after their last treatment in the pivotal study.

Subjects underwent repeat treatment with JUVEDERM® Ultra Plus in both NLFs. Demographics for subjects enrolled in the repeat treatment extended follow-up study were similar to those in the pivotal study. Routine follow-up visits for safety and effectiveness occurred at 4, 12, 24, 36, and 48 weeks after the repeat treatment. The investigator evaluated each subject for signs and symptoms of serious or unanticipated adverse events. The investigator also evaluated the severity of the subject's NLFs using the validated 5-point (range 0 to 4) photographic NLF severity scale. The subject made independent self-assessments of NLF severity using the non-photographic 5-point grading scale.

No serious or unanticipated adverse events were reported. The effectiveness results for repeat treatment with JUVEDERM® Ultra Plus based on the investigator's assessment of NLF severity after repeat treatment are presented in Table 8.

Table 8. Follow-up After Repeat Treatment Effectiveness Summary Investigator's NLF Severity Scores

	JUVEDERM® Ultra Plus (N = 24)	Improvement Since Baseline ^a
n ^b	24	2.7
Pre-repeat Treatment	24	2.2
Week 12	23	0.8
Week 24	23	1.0
Week 48	10	1.3
Improvement Since Baseline ^c	1.5	1.5

^a Number of subject NLFs with data at baseline and the specified time point

^b Mean score

Throughout the 48-week follow-up period, JUVEDERM® Ultra Plus provided a clinically significant improvement in NLF severity (≥ 1-point mean improvement) with a large majority of subjects treated with JUVEDERM® Ultra Plus demonstrating improvement at 24 weeks and beyond: 91% (21/23) at 24 weeks, and 90% (9/10) at 48 weeks (1 year).

C. Clinical Study for JUVEDERM® Ultra Plus XC

A prospective, double-blind, randomized, within-subject, controlled, multicenter clinical study was conducted to evaluate the safety and effectiveness of JUVEDERM® Ultra Plus XC compared with JUVEDERM® Ultra Plus without lidocaine. The purpose of this study was to evaluate the level of procedural pain (pain during injection) experienced by subjects when treated with each product. The duration of the study was 2 weeks.

A total of 36 subjects received a single treatment with JUVEDERM® Ultra Plus XC in one NLF and JUVEDERM® Ultra Plus without lidocaine in the other NLF. Within 30 minutes after both NLFs were treated, the subjects rated procedural pain on an 11-point scale and a 5-point comparative scale. Both the investigators and subjects rated NLF severity at baseline and 2 weeks after treatment using the 5-point NLF severity scale from the pivotal study. Subjects utilized an interactive voice-response-system diary to record common treatment site reactions for 14 days.

Most of the subjects were women (81% of Caucasian descent (75% with Fitzpatrick Skin Phototype II or III) (78%). Persons of color (Fitzpatrick Skin Phototypes IV, V, or VI) comprised 20% of treated subjects. Median age at study entry was 53 years (range, 32 to 80). Subject demographics are shown in Table 9.

Table 9. Subject Demographics (Number/% of Subjects)
N = 36 Subjects

Gender	Female	35	97%
Male	1	3%	
Ethnicity	29	81%	
Caucasian	5	14%	
African American	0	0%	
Hispanic	0	0%	
Other	1	3%	
Fitzpatrick Skin Type	I	1	3%
II	15	42%	
III	13	36%	
IV	2	6%	
V	2	6%	
VI	3	8%	

The pain scores for the NLFs treated with JUVEDERM® Ultra Plus XC were significantly lower (P < 0.0001) than for the NLFs treated with JUVEDERM® Ultra Plus without lidocaine (Table 10) based on the 11-point scale. On the comparative scale, 92% (33/36) of subjects rated the side with lidocaine as less or slightly less painful compared to the side with lidocaine as less or slightly less

Table 10. Subject Assessment of Procedural Pain Scores (N = 36)

	Mean Pain Score ^a
JUVEDERM® Ultra Plus XC	2.4
JUVEDERM® Ultra Plus	5.7
Mean Difference	-3.2

^a Procedural pain score ranges from 0 to 10 where 0 = No Pain and 10 = Worst Pain Imaginable

Table 11. Subject Assessments of Comparative Procedural Pain Score

	JUVEDERM® Ultra Plus (N (%))
JUVEDERM® Ultra Plus XC is less painful	28 (84%)
JUVEDERM® Ultra Plus XC is slightly less painful	10 (28%)
No difference between products	0 (0%)
JUVEDERM® Ultra Plus XC is slightly more painful	1 (3%)
JUVEDERM® Ultra Plus XC is more painful	2 (6%)

NLF severity improvement after 2 weeks was similar for both JUVEDERM® products (with and without lidocaine). Mean baseline score was 2.8, and a clinically significant improvement (severity reduction) to 1.0 was observed after 2 weeks for both products.


8. INSTRUCTIONS FOR USE

A. To Attach Needle to Syringe

STEP 1: Remove tip cap

Hold syringe and pull tip cap off the syringe as shown in Figure A.

Figure A



STEP 2: Insert needle

Hold the syringe body and firmly insert the hub of the needle (provided in the JUVEDERM® package) into the LUER-LOK® end of the syringe.

STEP 3: Tighten the needle

Tighten the needle by turning it firmly in a clockwise direction (see Figure B) until it is seated in the proper position as shown in Figure C.

Figure B




Figure C

Figure D

Figure X

Figure Y

Figure Z

Figure AA

Figure AB

Figure AC

Figure AD

Figure AE

Figure AF

Figure AG

Figure AH

Figure AI

Figure AJ

Figure AK

Figure AL

Figure AM

Figure AN

Figure AO

Figure AP

Figure AQ

Figure AR

Figure AS

Figure AT

Figure AU

Figure AV

Figure AW

Figure AX

Figure AY

Figure AZ

Figure BA

Figure BB

Figure BC

Figure BD

Figure BE

Figure BF

Figure BG

Figure BH

Figure BI

Figure BJ

Figure BK

Figure BL

Figure BM

Figure BN

Figure BO

Figure BP

Figure BQ

Figure BR

Figure BS

Figure BT

Figure BU

Figure BV

Figure BW

Figure BX

Figure BY

Figure BZ

Figure CA

Figure CB

Figure CC

Figure CD

Figure CE

Figure CF

Figure CG

Figure CH

Figure CI

Figure CJ

Figure CK

Figure CL

Figure CM

Figure CN

Figure CO

Figure CP

Figure CQ

Figure CR

Figure CS

Figure CT

Figure CU

Figure CV

Figure CW

Figure CX

Figure CY

Figure CZ

Figure DA

Figure DB

Figure DC

Figure DD

Figure DE

Figure DF

Figure DG

Figure DH

Figure DI

Figure DJ

Figure DK

Figure DL

Figure DM

Figure DN

Figure DO

Figure DP

Figure DQ

Figure DR

Figure DS

Figure DT

Figure DU

Figure DV

Figure DW

Figure DX

Figure DY

Figure DZ

Figure EA

Figure EB

Figure EC

Figure ED

Figure EE

Figure EF

Figure EG

Figure EH

Figure EI

Figure EJ

Figure EK

Figure EL

Figure EM

Figure EN

Figure EO

Figure EP

Figure EQ

Figure ER

Figure ES

Figure ET

Figure EU

Figure EV

Figure EW

Figure EX

Figure EY

Figure EZ

Figure FA

Figure FB

Figure FC

Figure FD

Figure FE

Figure FF

Figure FG

Figure FH

Figure FI

Figure FJ

Figure FK

Figure FL

Figure FM

Figure FN

Figure FO

Figure FP

Figure FQ

Figure FR


Figure FS

Figure FT

Figure FU

STEP 4: Remove the needle cap

Hold the syringe body in one hand and the needle cap in the other. Without twisting, pull in opposite directions to remove the needle cap as shown in Figure E.



B. Health Care Professional Instructions

1. JUVÉDERM® Ultra Plus XC injectable gel is a more highly cross-linked robust formulation, injected using a 27-G needle for volumizing and correction of deeper folds and wrinkles. Prior to treatment, the patient's medical history should be obtained and the patient should be fully apprised of the indications, contraindications, warnings, precautions, treatment responses, adverse reactions, and method of administration. Patients also should be advised that supplemental "touch-up" implantations may be required to achieve and maintain maximum correction.
2. The patient's soft-tissue deficiencies should be characterized with regard to etiology, distensibility, stress at the site, and depth of lesion. Depending on the type of skin, best results are obtained when the defect is readily distensible and correction can be visualized by manual manipulation (stretching) of the skin. Pretreatment photographs are recommended.
3. Although the study showed JUVÉDERM® Ultra Plus XC to be less painful than JUVÉDERM® Ultra Plus, supplementary anesthesia may be used for additional pain management during and after injection.
4. After ensuring that the patient has thoroughly washed the treatment area with soap and water, the area should be swabbed with alcohol or other antiseptic. Prior to injecting, depress the plunger rod until the product flows out of the needle.
5. After the first small amount of material has been injected into the patient, wait a full 3 seconds to allow the lidocaine to take effect before proceeding with the rest of the injection.
6. The injection technique may vary with regard to the angle and orientation of the bevel, the depth of injection, and the quantity administered. A linear threading technique, serial puncture injections, or a combination of the 2 have been used to achieve optimal results. Injecting the product too superficially may result in visible lumps and/or discoloration.
7. Inject JUVÉDERM® Ultra Plus XC applying even pressure on the plunger rod while slowly pulling the needle backward. The wrinkle should be lifted and eliminated by the end of the injection. It is important that the injection be stopped just before the needle is pulled out of the skin to prevent material from leaking out or ending up too superficially in the skin.
8. If the needle is blocked, do not increase the pressure on the plunger rod. Instead, stop the injection and replace the needle.
9. The typical total volume to achieve optimal correction of moderate to severe facial wrinkles and nasolabial folds is 1.6 mL per treatment site. The typical volume to achieve optimal correction for repeat treatment is 0.7 mL per treatment site.
10. Correct to 100% of the desired volume effect. Do not overcorrect. The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. Markedly indurated defects may be difficult to correct.
11. If immediate blanching occurs, the injection should be stopped and the area massaged until it returns to a normal color. Blanching may represent a vessel occlusion. If normal skin coloring does not return, do not continue with the injection.

Treat in accordance with American Society for Dermatologic Surgery guidelines, which include hyaluronidase injection.¹

12. When injection is completed, the treated site should be gently massaged so that it conforms to the contour of the surrounding tissues. If overcorrection occurs, massage the area between your fingers or against an underlying superficial bone to obtain optimal results.
13. With patients who have localized swelling, the degree of correction is sometimes difficult to judge at the time of treatment. In these cases, it is better to invite the patient to a touch-up session after 1-2 weeks.
14. Patients may have mild to moderate injection site responses, which typically resolve in a few days. If the treated area is swollen immediately after the injection, an ice pack can be applied to the site for a short period.
15. After the initial treatment, an additional treatment (from 1 to 2 weeks later) may be necessary to achieve the desired level of correction. If the wrinkle needs further treatment, the same procedure should be repeated until a satisfactory result is obtained. The need for an additional treatment may vary from patient to patient and is dependent upon a variety of factors such as wrinkle severity, skin elasticity, and dermal thickness at the treatment site.
16. The health care professional should instruct the patient to promptly report to her/him any evidence of problems possibly associated with the use of JUVÉDERM® Ultra Plus XC.

C. Patient Instructions

It is recommended that the following information be shared with patients:

- Within the first 24 hours, patients should avoid strenuous exercise, extensive sun or heat exposure, and alcoholic beverages. Exposure to any of the above may cause temporary redness, swelling, and/or itching at the injection sites.
- To report an adverse reaction, phone the Allergan® Product Support Department, 1-877-345-5372.

9. HOW SUPPLIED

JUVÉDERM® Ultra Plus XC injectable gel is supplied in individual treatment syringes with 27-G needles for single-patient use and ready for injection (implantation). The volume in each syringe is as stated on the syringe label and on the carton. The contents of the syringe are sterile and non-pyrogenic. Do not resterilize. Do not use if package is opened or damaged.

10. SHELF LIFE AND STORAGE

JUVÉDERM® Ultra Plus XC injectable gel must be used prior to the expiration date printed on the label.

Store at room temperature (up to 25°C/77°F). DO NOT FREEZE.

JUVÉDERM® Ultra Plus XC injectable gel has a clear appearance. In the event that a syringe contains material that is not clear, do not use the syringe; notify Allergan® Product Support immediately at 1-877-345-5372.

To place an order, contact Allergan® at 1-800-377-7790.

Allergan Aesthetics

an AbbVie company

Irvine, CA 92612 USA

1-800-624-4261

Made in France

© 2023 AbbVie. All rights reserved.

JUVÉDERM and its design are trademarks of Allergan Holdings France SAS, an AbbVie company, or its affiliates.

Patented. See: www.abbvie.com/patents

hcp.Juvederm.com

05/2023 022775 20079132

¹Alam M, Gladstone H, Kramer EM, et al. ASDS guidelines of care: injectable fillers. *Dermatol Surg.* 2008;34(suppl 1):S115-S148.